


|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

### Primary Tray

|    | Product #  | Product Desc.                                              | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|----|------------|------------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 1  | 04.118.510 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/10MM | 2   | 10886982108746 | <br>(01)10886982108746   |
| 2  | 04.118.512 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/12MM | 2   | 10886982108760 | <br>(01)10886982108760   |
| 3  | 04.118.514 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/14MM | 2   | 10886982108784 | <br>(01)10886982108784   |
| 4  | 04.118.516 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/16MM | 2   | 10886982108807 | <br>(01)10886982108807   |
| 5  | 04.118.518 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/18MM | 2   | 10886982108821 | <br>(01)10886982108821  |
| 6  | 04.118.520 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/20MM | 2   | 10886982108845 | <br>(01)10886982108845 |
| 7  | 04.118.522 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/22MM | 2   | 10886982108869 | <br>(01)10886982108869 |
| 8  | 04.118.524 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/24MM | 2   | 10886982108883 | <br>(01)10886982108883 |
| 9  | 04.118.526 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/26MM | 2   | 10886982108906 | <br>(01)10886982108906 |
| 10 | 04.118.528 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/28MM | 2   | 10886982108920 | <br>(01)10886982108920 |
| 11 | 04.118.530 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/30MM | 2   | 10886982108944 | <br>(01)10886982108944 |
| 12 | 04.118.532 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/32MM | 2   | 10886982108968 | <br>(01)10886982108968 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|    | Product #  | Product Desc.                                              | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|----|------------|------------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 13 | 04.118.534 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/34MM | 2   | 10886982108982 | <br>(01)10886982108982   |
| 14 | 04.118.536 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/36MM | 2   | 10886982109002 | <br>(01)10886982109002   |
| 15 | 04.118.538 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/38MM | 2   | 10886982109026 | <br>(01)10886982109026   |
| 16 | 04.118.540 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/40MM | 2   | 10886982109040 | <br>(01)10886982109040   |
| 17 | 04.118.542 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/42MM | 2   | 10886982109064 | <br>(01)10886982109064  |
| 18 | 04.118.544 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/44MM | 2   | 10886982109088 | <br>(01)10886982109088 |
| 19 | 04.118.546 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/46MM | 2   | 10886982109101 | <br>(01)10886982109101 |
| 20 | 04.118.548 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/48MM | 2   | 10886982109125 | <br>(01)10886982109125 |
| 21 | 04.118.550 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/50MM | 2   | 10886982109149 | <br>(01)10886982109149 |
| 22 | 04.118.552 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/52MM | 2   | 10886982109163 | <br>(01)10886982109163 |
| 23 | 04.118.554 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/54MM | 2   | 10886982109187 | <br>(01)10886982109187 |
| 24 | 04.118.556 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/56MM | 2   | 10886982109200 | <br>(01)10886982109200 |
| 25 | 04.118.558 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8 STRDRV RECESS/58MM | 2   | 10886982109224 | <br>(01)10886982109224 |


|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|    | Product #  | Product Desc.                                                   | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|----|------------|-----------------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 26 | 04.118.560 | 2.7MM TI METAPHYSEAL SCR/SLF- TPNG W/T8<br>STRDRV RECESS/60MM   | 2   | 10886982109248 | <br>(01)10886982109248   |
| 27 | 04.211.010 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 10MM | 3   | 10886982118950 | <br>(01)10886982118950   |
| 28 | 04.211.012 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 12MM | 3   | 10886982118974 | <br>(01)10886982118974   |
| 29 | 04.211.014 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 14MM | 3   | 10886982118998 | <br>(01)10886982118998   |
| 30 | 04.211.016 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 16MM | 3   | 10886982119001 | <br>(01)10886982119001  |
| 31 | 04.211.018 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 18MM | 3   | 10886982119018 | <br>(01)10886982119018 |
| 32 | 04.211.020 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 20MM | 3   | 10886982119025 | <br>(01)10886982119025 |
| 33 | 04.211.022 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 22MM | 3   | 10886982119032 | <br>(01)10886982119032 |
| 34 | 04.211.024 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 24MM | 3   | 10886982119049 | <br>(01)10886982119049 |
| 35 | 04.211.026 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 26MM | 3   | 10886982119056 | <br>(01)10886982119056 |
| 36 | 04.211.028 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 28MM | 3   | 10886982119063 | <br>(01)10886982119063 |
| 37 | 04.211.030 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 30MM | 3   | 10886982119070 | <br>(01)10886982119070 |
| 38 | 04.211.032 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8<br>STARDRIVE RECESS 32MM | 3   | 10886982119087 | <br>(01)10886982119087 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|    | Product #  | Product Desc.                                                | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|----|------------|--------------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 39 | 04.211.034 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 34MM | 3   | 10886982119094 | <br>(01)10886982119094   |
| 40 | 04.211.036 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 36MM | 3   | 10886982119100 | <br>(01)10886982119100   |
| 41 | 04.211.038 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 38MM | 3   | 10886982119117 | <br>(01)10886982119117   |
| 42 | 04.211.040 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 40MM | 3   | 10886982119124 | <br>(01)10886982119124   |
| 43 | 04.211.042 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 42MM | 3   | 10886982119131 | <br>(01)10886982119131  |
| 44 | 04.211.044 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 44MM | 3   | 10886982119148 | <br>(01)10886982119148 |
| 45 | 04.211.046 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 46MM | 3   | 10886982119155 | <br>(01)10886982119155 |
| 46 | 04.211.048 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 48MM | 3   | 10886982119162 | <br>(01)10886982119162 |
| 47 | 04.211.050 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 50MM | 3   | 10886982119179 | <br>(01)10886982119179 |
| 48 | 04.211.052 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 52MM | 3   | 10886982119186 | <br>(01)10886982119186 |
| 49 | 04.211.054 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 54MM | 3   | 10886982119193 | <br>(01)10886982119193 |
| 50 | 04.211.056 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 56MM | 3   | 10886982119209 | <br>(01)10886982119209 |
| 51 | 04.211.058 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 58MM | 3   | 10886982119216 | <br>(01)10886982119216 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|    | Product #  | Product Desc.                                                | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|----|------------|--------------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 52 | 04.211.060 | 2.7MM TI VA LCKNG SCR SLF-TPNG WITH T8 STARDRIVE RECESS 60MM | 3   | 10886982119223 | <br>(01)10886982119223   |
| 53 | 04.212.054 | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 54MM     | 3   | 10886982120618 | <br>(01)10886982120618   |
| 54 | 04.212.056 | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 56MM     | 3   | 10886982120625 | <br>(01)10886982120625   |
| 55 | 04.212.058 | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 58MM     | 3   | 10886982120632 | <br>(01)10886982120632   |
| 56 | 402.210    | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-10MM  | 3   | 10886982205728 | <br>(01)10886982205728  |
| 57 | 402.212    | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-12MM  | 3   | 10886982205742 | <br>(01)10886982205742 |
| 58 | 402.214    | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-14MM  | 3   | 10886982205766 | <br>(01)10886982205766 |
| 59 | 402.216    | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-16MM  | 3   | 10886982205773 | <br>(01)10886982205773 |
| 60 | 402.218    | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-18MM  | 3   | 10886982205780 | <br>(01)10886982205780 |
| 61 | 402.220    | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-20MM  | 3   | 10886982205797 | <br>(01)10886982205797 |
| 62 | 402.222    | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-22MM  | 3   | 10886982205803 | <br>(01)10886982205803 |
| 63 | 402.224    | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-24MM  | 3   | 10886982205810 | <br>(01)10886982205810 |
| 64 | 402.226    | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-26MM  | 3   | 10886982205827 | <br>(01)10886982205827 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|    | Product # | Product Desc.                                               | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|----|-----------|-------------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 65 | 402.228   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-28MM | 3   | 10886982205834 | <br>(01)10886982205834   |
| 66 | 402.230   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-30MM | 3   | 10886982205841 | <br>(01)10886982205841   |
| 67 | 402.232   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-32MM | 3   | 10886982205858 | <br>(01)10886982205858   |
| 68 | 402.234   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-34MM | 3   | 10886982205865 | <br>(01)10886982205865   |
| 69 | 402.236   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-36MM | 3   | 10886982205872 | <br>(01)10886982205872  |
| 70 | 402.238   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-38MM | 3   | 10886982205889 | <br>(01)10886982205889 |
| 71 | 402.240   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-40MM | 3   | 10886982205896 | <br>(01)10886982205896 |
| 72 | 402.242   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-42MM | 3   | 10886982205902 | <br>(01)10886982205902 |
| 73 | 402.244   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-44MM | 3   | 10886982205919 | <br>(01)10886982205919 |
| 74 | 402.246   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-46MM | 3   | 10886982205926 | <br>(01)10886982205926 |
| 75 | 402.248   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-48MM | 3   | 10886982205933 | <br>(01)10886982205933 |
| 76 | 402.250   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-50MM | 3   | 10886982205940 | <br>(01)10886982205940 |
| 77 | 402.255   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-55MM | 3   | 10886982205957 | <br>(01)10886982205957 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|    | Product # | Product Desc.                                                | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|----|-----------|--------------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 78 | 402.260   | 2.7MM TI LOCKING SCR SLF-TPNG WITH T8 STARDRIVE RECESS-60MM  | 3   | 10886982205964 | <br>(01)10886982205964   |
| 79 | 402.870   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 10MM | 2   | 10886982206695 | <br>(01)10886982206695   |
| 80 | 402.872   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 12MM | 2   | 10886982206701 | <br>(01)10886982206701   |
| 81 | 402.874   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 14MM | 2   | 10886982206718 | <br>(01)10886982206718   |
| 82 | 402.876   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 16MM | 2   | 10886982206725 | <br>(01)10886982206725  |
| 83 | 402.878   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 18MM | 2   | 10886982206732 | <br>(01)10886982206732 |
| 84 | 402.880   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 20MM | 2   | 10886982206749 | <br>(01)10886982206749 |
| 85 | 402.882   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 22MM | 2   | 10886982206756 | <br>(01)10886982206756 |
| 86 | 402.884   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 24MM | 2   | 10886982206763 | <br>(01)10886982206763 |
| 87 | 402.886   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 26MM | 2   | 10886982206770 | <br>(01)10886982206770 |
| 88 | 402.888   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 28MM | 2   | 10886982206787 | <br>(01)10886982206787 |
| 89 | 402.890   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 30MM | 2   | 10886982206794 | <br>(01)10886982206794 |
| 90 | 402.892   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 32MM | 2   | 10886982206800 | <br>(01)10886982206800 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|     | Product # | Product Desc.                                                | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|-----|-----------|--------------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 91  | 402.894   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 34MM | 2   | 10886982206817 | <br>(01)10886982206817   |
| 92  | 402.896   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 36MM | 2   | 10886982206824 | <br>(01)10886982206824   |
| 93  | 402.898   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 38MM | 2   | 10886982206831 | <br>(01)10886982206831   |
| 94  | 402.900   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 40MM | 2   | 10886982206848 | <br>(01)10886982206848   |
| 95  | 402.962   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 42MM | 2   | 10886982206855 | <br>(01)10886982206855  |
| 96  | 402.963   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 44MM | 2   | 10886982206862 | <br>(01)10886982206862 |
| 97  | 402.965   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 46MM | 2   | 10886982206879 | <br>(01)10886982206879 |
| 98  | 402.966   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 48MM | 2   | 10886982206886 | <br>(01)10886982206886 |
| 99  | 402.967   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 50MM | 2   | 10886982206893 | <br>(01)10886982206893 |
| 100 | 402.968   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 55MM | 2   | 10886982206909 | <br>(01)10886982206909 |
| 101 | 402.969   | 2.7MM TI CORTEX SCREW SLF-TPNG WITH T8 STARDRIVE RECESS 60MM | 2   | 10886982206916 | <br>(01)10886982206916 |
| 102 | 404.810   | 3.5MM TI CORTEX SCREW SELF-TAPPING 10MM                      | 6   | 10886982207272 | <br>(01)10886982207272 |
| 103 | 404.812   | 3.5MM TI CORTEX SCREW SELF-TAPPING 12MM                      | 6   | 10886982207289 | <br>(01)10886982207289 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|     | Product # | Product Desc.                           | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|-----|-----------|-----------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 104 | 404.814   | 3.5MM TI CORTEX SCREW SELF-TAPPING 14MM | 6   | 10886982207296 | <br>(01)10886982207296   |
| 105 | 404.816   | 3.5MM TI CORTEX SCREW SELF-TAPPING 16MM | 6   | 10886982207302 | <br>(01)10886982207302   |
| 106 | 404.818   | 3.5MM TI CORTEX SCREW SELF-TAPPING 18MM | 6   | 10886982207319 | <br>(01)10886982207319   |
| 107 | 404.820   | 3.5MM TI CORTEX SCREW SELF-TAPPING 20MM | 6   | 10886982207326 | <br>(01)10886982207326   |
| 108 | 404.822   | 3.5MM TI CORTEX SCREW SELF-TAPPING 22MM | 6   | 10886982207333 | <br>(01)10886982207333  |
| 109 | 404.824   | 3.5MM TI CORTEX SCREW SELF-TAPPING 24MM | 6   | 10886982207340 | <br>(01)10886982207340 |
| 110 | 404.826   | 3.5MM TI CORTEX SCREW SELF-TAPPING 26MM | 6   | 10886982207357 | <br>(01)10886982207357 |
| 111 | 404.828   | 3.5MM TI CORTEX SCREW SELF-TAPPING 28MM | 6   | 10886982207364 | <br>(01)10886982207364 |
| 112 | 404.830   | 3.5MM TI CORTEX SCREW SELF-TAPPING 30MM | 6   | 10886982207371 | <br>(01)10886982207371 |
| 113 | 404.832   | 3.5MM TI CORTEX SCREW SELF-TAPPING 32MM | 3   | 10886982207388 | <br>(01)10886982207388 |
| 114 | 404.834   | 3.5MM TI CORTEX SCREW SELF-TAPPING 34MM | 3   | 10886982207395 | <br>(01)10886982207395 |
| 115 | 404.836   | 3.5MM TI CORTEX SCREW SELF-TAPPING 36MM | 3   | 10886982207401 | <br>(01)10886982207401 |
| 116 | 404.838   | 3.5MM TI CORTEX SCREW SELF-TAPPING 38MM | 3   | 10886982207418 | <br>(01)10886982207418 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|     | Product # | Product Desc.                                      | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|-----|-----------|----------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 117 | 404.840   | 3.5MM TI CORTEX SCREW SELF-TAPPING 40MM            | 3   | 10886982207425 | <br>(01)10886982207425   |
| 118 | 404.845   | 3.5MM TI CORTEX SCREW SELF-TAPPING 45MM            | 3   | 10886982207432 | <br>(01)10886982207432   |
| 119 | 404.850   | 3.5MM TI CORTEX SCREW SELF-TAPPING 50MM            | 3   | 10886982207449 | <br>(01)10886982207449   |
| 120 | 404.855   | 3.5MM TI CORTEX SCREW SELF-TAPPING 55MM            | 3   | 10886982207456 | <br>(01)10886982207456   |
| 121 | 406.010   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/10MM | 2   | 10886982207463 | <br>(01)10886982207463  |
| 122 | 406.012   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/12MM | 2   | 10886982207470 | <br>(01)10886982207470 |
| 123 | 406.014   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/14MM | 2   | 10886982207487 | <br>(01)10886982207487 |
| 124 | 406.016   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/16MM | 2   | 10886982207494 | <br>(01)10886982207494 |
| 125 | 406.018   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/18MM | 2   | 10886982207500 | <br>(01)10886982207500 |
| 126 | 406.020   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/20MM | 2   | 10886982207517 | <br>(01)10886982207517 |
| 127 | 406.022   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/22MM | 2   | 10886982207524 | <br>(01)10886982207524 |
| 128 | 406.024   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/24MM | 2   | 10886982207531 | <br>(01)10886982207531 |
| 129 | 406.026   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/26MM | 2   | 10886982207548 | <br>(01)10886982207548 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|     | Product # | Product Desc.                                            | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|-----|-----------|----------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 130 | 406.028   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/28MM       | 2   | 10886982207555 | <br>(01)10886982207555   |
| 131 | 406.030   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/30MM       | 2   | 10886982207562 | <br>(01)10886982207562   |
| 132 | 406.035   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/35MM       | 2   | 10886982207579 | <br>(01)10886982207579   |
| 133 | 406.040   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/40MM       | 2   | 10886982207586 | <br>(01)10886982207586   |
| 134 | 406.045   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/45MM       | 2   | 10886982207593 | <br>(01)10886982207593  |
| 135 | 406.050   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/50MM       | 2   | 10886982207609 | <br>(01)10886982207609 |
| 136 | 406.055   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/55MM       | 2   | 10886982207616 | <br>(01)10886982207616 |
| 137 | 406.060   | 4.0MM TI CANCELLOUS BONE SCREW FULLY THREADED/60MM       | 2   | 10886982207623 | <br>(01)10886982207623 |
| 138 | 412.101   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 10MM | 3   | 10886982210340 | <br>(01)10886982210340 |
| 139 | 412.102   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 12MM | 3   | 10886982210357 | <br>(01)10886982210357 |
| 140 | 412.103   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 14MM | 3   | 10886982210364 | <br>(01)10886982210364 |
| 141 | 412.104   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 16MM | 3   | 10886982210371 | <br>(01)10886982210371 |
| 142 | 412.105   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 18MM | 3   | 10886982210388 | <br>(01)10886982210388 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|     | Product # | Product Desc.                                            | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|-----|-----------|----------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 143 | 412.106   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 20MM | 3   | 10886982210395 | <br>(01)10886982210395   |
| 144 | 412.107   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 22MM | 3   | 10886982210401 | <br>(01)10886982210401   |
| 145 | 412.108   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 24MM | 3   | 10886982210418 | <br>(01)10886982210418   |
| 146 | 412.109   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 26MM | 3   | 10886982210425 | <br>(01)10886982210425   |
| 147 | 412.110   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 28MM | 3   | 10886982210449 | <br>(01)10886982210449  |
| 148 | 412.111   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 30MM | 3   | 10886982210463 | <br>(01)10886982210463 |
| 149 | 412.112   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 32MM | 3   | 10886982210487 | <br>(01)10886982210487 |
| 150 | 412.113   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 34MM | 3   | 10886982210500 | <br>(01)10886982210500 |
| 151 | 412.115   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 36MM | 3   | 10886982210548 | <br>(01)10886982210548 |
| 152 | 412.116   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 38MM | 3   | 10886982210562 | <br>(01)10886982210562 |
| 153 | 412.117   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 40MM | 3   | 10886982210586 | <br>(01)10886982210586 |
| 154 | 412.118   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 42MM | 3   | 10886982210609 | <br>(01)10886982210609 |
| 155 | 412.119   | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 45MM | 3   | 10886982210623 | <br>(01)10886982210623 |

|                 |                                                                                                         |
|-----------------|---------------------------------------------------------------------------------------------------------|
| Set Code:       | 01.133.408E                                                                                             |
| Set Barcode:    | <br>(01)10886982287120 |
| Description:    | TI - UNIVERSAL SMALL FRAGMENT SCREW SET EVAL                                                            |
| Case Report ID: |                                                                                                         |

 **Primary Tray**

|     | Product #  | Product Desc.                                               | QTY | UDI-DI         | UDI Barcode                                                                                                 |
|-----|------------|-------------------------------------------------------------|-----|----------------|-------------------------------------------------------------------------------------------------------------|
| 156 | 412.120    | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 48MM    | 3   | 10886982210647 | <br>(01)10886982210647   |
| 157 | 412.121    | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 50MM    | 3   | 10886982210661 | <br>(01)10886982210661   |
| 158 | 412.122    | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 52MM    | 3   | 10886982210685 | <br>(01)10886982210685   |
| 159 | 412.123    | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 55MM    | 3   | 10886982210708 | <br>(01)10886982210708   |
| 160 | 412.124    | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 60MM    | 3   | 10886982210722 | <br>(01)10886982210722  |
| 161 | 412.134    | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 44MM    | 3   | 10886982210814 | <br>(01)10886982210814 |
| 162 | 412.136    | 3.5MM TI LOCKING SCR SLF-TPNG WITH STARDRIVE RECESS 46MM    | 3   | 10886982210838 | <br>(01)10886982210838 |
| 163 | 419.98     | WASHER 7.0MM                                                | 10  | 10886982214829 | <br>(01)10886982214829 |
| 164 | 60.116.507 | SCREW TYPE PUSH PIN/CORTEX FOR MODULAR GRAPHIC CASE SYS     | 1   | 20886982246551 | <br>(01)20886982246551 |
| 165 | 60.116.513 | SCREW TYPE PUSH PIN/LOCKING FOR MODULAR GRAPHIC CASE SYS    | 1   | 20886982246612 | <br>(01)20886982246612 |
| 166 | 60.116.521 | SCREW TYPE PUSH PIN/VA LOCKING FOR MODULAR GRAPHIC CASE SYS | 1   | 20886982246698 | <br>(01)20886982246698 |
| 167 | 60.116.527 | SCR TYPE PUSH PIN/METAPHYSEAL FOR MODULAR GRAPHIC CASE SYS  | 2   | 10886982246752 | <br>(01)10886982246752 |
| 168 | 60.133.150 | UNIVERSAL SF SCREW RACK                                     | 1   | 10886982269751 | <br>(01)10886982269751 |