



Set Code:	115.720ACL
Set Barcode:	 (01)10886982311733
Description:	LARGE EXTERNAL FIXATOR SET W/ SELF-DRILLING SCHANZ SCREWS
Case Report ID:	

 Primary Tray

	Product #	Product Desc.	QTY	UDI-DI	UDI Barcode
1	293.74	5.0MM STEINMANN PIN WITH CENTRAL THREAD 200MM	4	10886982181282	 (01)10886982181282
2	294.56	5.0MM SCHANZ SCREW BLUNTED TROCAR POINT 200MM	8	10886982181572	 (01)10886982181572
3	294.784	5.0MM SELF-DRILLING SCHANZ SCREW 60MM THRD/150MM	4	10886982182135	 (01)10886982182135
4	294.785	5.0MM SELF-DRILLING SCHANZ SCREW 60MM THRD/175MM	8	10886982182166	 (01)10886982182166
5	294.786	5.0MM SELF-DRILLING SCHANZ SCREW 80MM THRD/200MM	8	10886982182197	 (01)10886982182197
6	294.950	6.0MM TRANSFIXATION PIN 225MM	4	10886982182432	 (01)10886982182432
7	310.37	3.5MM DRILL BIT/QC/195MM	2	10886982186553	 (01)10886982186553
8	310.48	4.5MM DRILL BIT/QC/195MM	2	10886982186683	 (01)10886982186683
9	321.20	RATCHET WRENCH-11MM WIDTH ACROSS FLATS	2	10886982190444	 (01)10886982190444
10	355.14	CANNULATED SOCKET WRENCH- 11MM WIDTH ACROSS FLATS	1	10886982194800	 (01)10886982194800
11	390.002	LARGE EX-FIX MULTI-PIN CLAMP MR-CONDITIONAL / 6-POSITION	4	10886982197795	 (01)10886982197795
12	390.003	LARGE EX-FIX ROD ATTCHMT/MULTI PIN CLAMP/MR-CONDITIONAL	6	10886982197801	 (01)10886982197801

Set Code:	115.720ACL
Set Barcode:	 (01)10886982311733
Description:	LARGE EXTERNAL FIXATOR SET W/ SELF-DRILLING SCHANZ SCREWS
Case Report ID:	

 Primary Tray

	Product #	Product Desc.	QTY	UDI-DI	UDI Barcode
13	390.004	LARGE EX-FIX MULTI-PIN CLAMP MR-CONDITIONAL / 4-POSITION	2	10886982197818	 (01)10886982197818
14	390.005	LARGE EX-FIX COMBINATION CLAMP MR-CONDITIONAL	12	10886982197825	 (01)10886982197825
15	390.006	LRG EX-FIX DYNAMIZATION CLAMP F/COMBO CLAMP / MR-CONDITIONAL	4	10886982197832	 (01)10886982197832
16	390.007	LARGE EX-FIX TUBE-TUBE CLAMP MR-CONDITIONAL	2	10886982197849	 (01)10886982197849
17	390.008	LARGE EX-FIX OPEN ADJ CLAMP MR-CONDITIONAL	8	10886982197856	 (01)10886982197856
18	392.951	8.0MM/6.0MM THREADED DRILL SLEEVE-SHORT	1	10886982198822	 (01)10886982198822
19	392.952	8.0MM/6.0MM THREADED DRILL SLEEVE-LONG	1	10886982198839	 (01)10886982198839
20	392.963	6-POSITION DRILL GUIDE HANDLE	1	10886982198884	 (01)10886982198884
21	393.10	UNIVERSAL CHUCK WITH T-HANDLE	1	10886982198983	 (01)10886982198983
22	393.103	DRIVE ADAPTOR WITH QC FOR 5.0MM SCHANZ SCREWS	1	10886982199010	 (01)10886982199010
23	393.104	DRIVE ADAPTOR WITH QC FOR 6.0MM SCHANZ SCREWS	1	10886982199027	 (01)10886982199027
24	393.746	SPLIT TISSUE PROTECTION SLEEVE 5.0MM	1	10886982199393	 (01)10886982199393
25	393.76	OPEN COMPRESSOR	2	10886982199430	 (01)10886982199430

Set Code:	115.720ACL
Set Barcode:	 (01)10886982311733
Description:	LARGE EXTERNAL FIXATOR SET W/ SELF-DRILLING SCHANZ SCREWS
Case Report ID:	

 Primary Tray

	Product #	Product Desc.	QTY	UDI-DI	UDI Barcode
26	394.181	3.5MM TROCAR-SHORT	1	10886982199621	 (01)10886982199621
27	394.182	3.5MM TROCAR-LONG	1	10886982199638	 (01)10886982199638
28	394.80	LARGE EX-FIX 11MM CRBN FBR ROD 100MM / MR-CONDITIONAL	4	10886982199850	 (01)10886982199850
29	394.82	LARGE EX-FIX 11MM CRBN FBR ROD 150MM / MR-CONDITIONAL	4	10886982199874	 (01)10886982199874
30	394.83	LARGE EX-FIX 11MM CRBN FBR ROD 200MM / MR-CONDITIONAL	4	10886982199881	 (01)10886982199881
31	394.84	LARGE EX-FIX 11MM CRBN FBR ROD 250MM / MR-CONDITIONAL	4	10886982199898	 (01)10886982199898
32	394.85	LARGE EX-FIX 11MM CRBN FBR ROD 300MM / MR-CONDITIONAL	4	10886982199904	 (01)10886982199904
33	394.86	LARGE EX-FIX 11MM CRBN FBR ROD 350MM / MR-CONDITIONAL	4	10886982199911	 (01)10886982199911
34	394.87	LARGE EX-FIX 11MM CRBN FBR ROD 400MM / MR-CONDITIONAL	4	10886982199928	 (01)10886982199928
35	394.97	PROTECTIVE CAPS FOR 11.0MM TUBES & CARBON FIBER RODS	1	20886982199987	 (01)20886982199987
36	394.993	PROTECTIVE CAPS FOR 5.0MM FIXATION PINS	1	20886982200065	 (01)20886982200065
37	394.994	PROTECTIVE CAPS FOR 6.0MM FIXATION PINS	1	20886982200089	 (01)20886982200089
38	395.911	DRILL SLEEVE HANDLE	1	10886982200891	 (01)10886982200891

Set Code:	115.720ACL
Set Barcode:	 (01)10886982311733
Description:	LARGE EXTERNAL FIXATOR SET W/ SELF-DRILLING SCHANZ SCREWS
Case Report ID:	

 Primary Tray

	Product #	Product Desc.	QTY	UDI-DI	UDI Barcode
39	395.912	5.0MM/3.5MM DRILL SLEEVE-SHORT	1	10886982200907	 (01)10886982200907
40	395.913	5.0MM/3.5MM DRILL SLEEVE-LONG	1	10886982200914	 (01)10886982200914
41	395.921	6.0MM/5.0MM THREADED DRILL SLEEVE-SHORT	1	10886982200921	 (01)10886982200921
42	395.923	6.0MM/5.0MM THREADED DRILL SLEEVE-LONG	1	10886982200945	 (01)10886982200945
43	690.315	LARGE EXTERNAL FIXATOR GRAPHIC CASE	1	10886982252104	 (01)10886982252104